

Job Aid: Ultrasound Pelvis-Non OB Protocol		
Division: SLH,SNGH	Manual:	Locations: Radiology Ultrasound
Original Date: 8/1/12	Section:	
Owner: Imaging Operations		
Approved By: Richard Thomas, MD, Department Chief		Revision Date: 11/20/13, 3/04/14, 1/20/17, 2/06/2020

TRANSABDOMINAL PELVIC ULTRASOUND

Organs demonstrated: Uterus, endometrium, adnexae, and ovaries

- | | |
|---|---|
| <p>1. <u>Sagittal Uterus</u>
 1 image midline w/wo calipers
 1 image midline w/wo calipers
 3 images (R1-R3)
 3 images (L1-L3)</p> | <p>3. <u>RT Adnexa/Ovary</u>
 sagittal ovary w/wo calipers
 TRV ovary w/wo calipers
 color and pulsed Doppler
 demonstrating perfusion</p> |
| <p>2. <u>Transverse Uterus</u>
 1 image vagina
 1 image cervix
 1 image uterus body
 1 image uterine fundus w/wo calipers</p> | <p>4. <u>LT Adnexa/Ovary</u>
 sagittal ovary w/wo calipers
 TRV ovary w/wo calipers
 color and pulsed Doppler
 demonstrating perfusion</p> |

Imaging Notes:

Do not need to measure all fibroids...Measure largest 1-3 fibroids or any designated for treatment. Show any submucosal and pedunculated fibroids.

When measuring fibroids and/or masses use Dual screen or make sure Sag and Trv Images are consecutive.

Measure Uterus from cervix to fundus...If prescription requests uterine volume- then measurement should be from internal os to fundus (do not include cervix)!

Related Documents:

<i>Policy</i>	N/A
<i>Procedure</i>	N/A